




ATTACHMENT A


Exemplary Infringement of U.S. Patent No. 10,227,639 by NanoString's CosMx Spatial Molecular Imaging Platform and Workflow


Claim No.	Claim Language	NanoString Spatial Molecular Imaging
1 (Preamble)	A method for analyte identification, comprising:	<p>To the extent that the preamble is found to be limiting, NanoString's CosMx Spatial Molecular Imaging ("CosMx SMI") workflow, including when used by NanoString through its CosMx Technology Access Program, is a method for analyte (RNA or protein) identification.</p> <p style="padding-left: 40px;">"SMI is a completely automated and integrated platform comprising chemistry, hardware, and software that enables highly sensitive spatial profiling of 980 RNAs and 108 proteins in FFPE tissues at single-cell and subcellular resolution."</p> <p style="padding-left: 40px;">"High-Plex Multiomic Analysis in FFPE at Single-Cellular and Subcellular Level by Spatial Molecular Imaging," Shanshan He <i>et al.</i>, (preprint available at https://www.biorxiv.org/content/10.1101/2021.11.03.467020v2) ("He-Beechem") at 4</p> <p style="padding-left: 40px;">"CosMx SMI is an integrated system with mature cyclic fluorescent in situ hybridization (FISH) chemistry, high-resolution imaging readout, and interactive data analysis and visualization software."</p>

Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		<div data-bbox="800 235 1890 560">  <p>1 SAMPLE PREPARATION</p> <p>Permeabilize, fix, retrieve targets Hybridization of RNA specific probes and antibodies Flow cell assembly</p> <p>Compatible with formalin-fixed paraffin embedded (FFPE) and fresh frozen (FF) tissues</p> <p>2 INTEGRATED READOUT</p> <p>CosMx Spatial Molecular Imager</p> <p>Robust <i>in situ</i> hybridization chemistry and readout</p> <p>3 INTERACTIVE DATA ANALYSIS</p> <p>Cloud-based scalable computing and storage with interactive data viewer</p> </div> <div data-bbox="814 581 1690 620"> <p>Easy Sample Preparation, Compatible with Any Sample Type</p> </div> <div data-bbox="842 649 1837 1039">  <p>START STEP 1 STEP 2 STEP 3</p> <p>FFPE tissue on standard slides Permeabilize, fix — Retrieve targets Hybridization- RNA specific probes and antibodies bind to target Assemble into flow cell</p> <p>Streamlined and simple workflow that integrates with standard ISH protocol with no need for tissue expansion or clearing, cDNA synthesis or amplification. Go from sample to result faster.</p> </div>


Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		<p data-bbox="842 240 1545 272">Automated Cyclic <i>in situ</i> Hybridization Chemistry</p>  <p data-bbox="852 638 1709 688">Robust hybridization chemistry that provides higher sensitivity and supports high plex assays in your tissue samples to uncover deeper biological insights.</p> <p data-bbox="1077 748 1896 857">“CosMx Spatial Molecular Imager,” available at https://www.nanostring.com/products/cosmx-spatial-molecular-imager/overview/ (“CosMx SMI Overview”)</p> <p data-bbox="884 894 1896 1073">“NanoString provides complete spatial molecular imaging service through Technology Access Program (TAP) from our state-of-the-art laboratory. The Spatial Molecular Imaging TAP enables researchers to visualize and quantify up to 1,000 RNA expression directly in single cells with their spatial context preserved.”</p> <p data-bbox="1077 1117 1896 1226">“CosMx SMI Technology Access Program,” https://www.nanostring.com/products/cosmx-spatial-molecular-imager/technology-access-program/ (“CosMx SMI TAP”)</p>

Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		<p>“How It Works:</p> <p>01: Submit your Formalin-Fixed Paraffin-Embedded tissue samples (whole mount tissue sections or tissue microarrays) to NanoString’s Translational Services Lab</p> <p>02: Nanostring will analyze your samples using spatial molecular imaging system with up to 1,000-plex RNA panel.</p> <p>03: Receive detailed report including raw data and analyzed results in 6-10 weeks.”</p> <p style="text-align: right;">CosMx SMI TAP</p>
1[a]	<p>(a) contacting a sample with a plurality of detection reagents, wherein said plurality of detection reagents comprises a detection reagent that targets an analyte of a plurality of analytes immobilized in the sample, wherein said detection reagent comprises:</p> <p>(i) probe targeting said analyte and</p> <p>(ii) a nucleic acid label comprising a plurality of pre-determined subsequences,</p>	<p>NanoString’s CosMx SMI workflow, including when performed by NanoString as part of its CosMx Technology Access Program, includes contacting a sample with a plurality of detection reagents, where the plurality of detection reagents comprises: a detection reagent that targets an analyte from among a plurality of analytes immobilized in the sample and includes a probe (RNA-specific probe or antibody probe) that targets the analyte (target RNA or protein, respectively) and that is conjugated together with a nucleic acid label comprising two or more of pre-determined subsequences comprising an identifier of the probe.</p>


Claim No.	Claim Language	NanoString Spatial Molecular Imaging
	<p>wherein said probe and said nucleic acid label are conjugated together, and wherein said plurality of pre-determined subsequences forms an identifier of said probe;</p>	<p>Easy Sample Preparation, Compatible with Any Sample Type</p>  <p>Streamlined and simple workflow that integrates with standard ISH protocol with no need for tissue expansion or clearing, cDNA synthesis or amplification. Go from sample to result faster.</p> <p>CosMx SMI Overview</p> <p>“For sample preparation, SMI utilizes standard methods performed for FISH on FFPE tissue sections to expose RNA targets, followed by the introduction of fluorescent bead-based fiducials that are fixed to the tissue to provide an optical reference for cyclic image acquisition and registration.”</p> <p>He-Beechem at 5</p>


Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		 <p data-bbox="1459 925 1896 958">He-Beechem, Fig. 1A (annotated)</p> <p data-bbox="884 998 1896 1250">“Schematic description of SMI ISH probe and reporter design. ISH probe consists of the target-binding domain and readout domain. The target-binding domain is a 35-50 nt DNA sequence that hybridizes with target RNA. The readout domain is a 60-80 nt DNA sequence that contains 4 consecutive 10-20 nt reporter-landing domains, where each landing domain can be hybridized with a unique reporter. With a 64-bit barcoding design, there are a total of 64 unique reporter-landing sequences.”</p> <p data-bbox="1220 1299 1896 1331">He-Beechem at 8, notes at Fig. 1A (emphasis added)</p>


Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		<p>He-Beechem, Fig. S13 (annotated)</p> <p>“Protein Detection Probe. The protein detection probe is readily adapted from the RNA detection probe shown in Figure 1. For protein detection, the readout domain contains an additional universal handle and is conjugated to the antibody via a site-specific linker.”</p> <p>He-Beechem at 53, notes at Fig. S13</p>

Claim No.	Claim Language	NanoString Spatial Molecular Imaging
1[b]	<p>(b) with said analyte immobilized in the sample and said probe coupled to said analyte,</p> <p>(i) hybridizing a first decoder probe with a first subsequence of said plurality of pre-determined subsequences, wherein said first decoder probe comprises a first detectable label,</p> <p>(ii) detecting a first signal signature from said first detectable label,</p> <p>(iii) hybridizing a second decoder probe with a second subsequence of said plurality of pre-determined subsequences, wherein said second decoder probe comprises a second detectable label, and</p> <p>(iv) detecting a second signal signature from said second detectable label, to provide a set of signal signatures comprising said</p>	<p>In NanoString's CosMx SMI workflow, including when performed by NanoString as part of its CosMx Technology Access Program, the analyte (RNA and/or protein) is immobilized in the sample and the probe (RNA-specific probe or antibody probe) is coupled (hybridized) to the analyte. A first decoder probe (reporter probe) comprising a first detectable (color) label is hybridized with a first subsequence of the label and a first signal signature is detected from the first detectable label. A second decoder probe comprising a second detectable label is hybridized with a second subsequence of the label and a second signal signature is detected from the second detectable label. A set of signal signatures is provided comprising the first signal signature and second signal signature.</p> <div data-bbox="835 683 1705 719"> <p>Easy Sample Preparation, Compatible with Any Sample Type</p>  <p>Streamlined and simple workflow that integrates with standard ISH protocol with no need for tissue expansion or clearing, cDNA synthesis or amplification. Go from sample to result faster.</p> <p style="text-align: right;">CosMx SMI Overview</p> </div> <p>“For sample preparation, SMI utilizes standard methods performed for FISH on FFPE tissue sections to expose RNA targets, followed by the introduction of</p>


Claim No.	Claim Language	NanoString Spatial Molecular Imaging
	first signal signature and said second signal signature; and	<p data-bbox="884 237 1898 305">fluorescent bead-based fiducials that are fixed to the tissue to provide an optical reference for cyclic image acquisition and registration.”</p> <p data-bbox="1671 354 1898 386">He-Beechem at 5</p> <div data-bbox="810 467 1843 938"> <p data-bbox="810 521 1031 662">a first decoder probe (reporter) with detectable (color) label</p> <p data-bbox="1031 472 1304 521">Reporter PC-Cleavable sites</p> <p data-bbox="1346 537 1465 565">ISH probe</p> <p data-bbox="1346 618 1444 662">Readout domain 60-80 nt</p> <p data-bbox="1535 699 1843 808">a second subsequence (second decoder probe not pictured)</p> <p data-bbox="1094 748 1262 792">Target binding domain 35-50 nt</p> <p data-bbox="1052 829 1108 857">RNA</p> <p data-bbox="1440 862 1612 938">a first subsequence</p> </div> <p data-bbox="1461 1003 1898 1036">He-Beechem, Fig. 1A (annotated)</p> <p data-bbox="884 1073 1898 1328">“Schematic description of SMI ISH probe and reporter design. ISH probe consists of the target-binding domain and readout domain. The target-binding domain is a 35-50 nt DNA sequence that hybridizes with target RNA. The readout domain is a 60-80 nt DNA sequence that contains 4 consecutive 10-20 nt reporter-landing domains, where each landing domain can be hybridized with a unique reporter. With a 64-bit barcoding design, there are a total of 64 unique reporter-landing sequences.”</p> <p data-bbox="1220 1382 1898 1414">He-Beechem at 8, notes at Fig. 1A (emphasis added)</p>

Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		 <p data-bbox="840 389 1050 535">a first decoder probe (reporter) with detectable (color) label</p> <p data-bbox="1029 341 1386 422">Reporter PC-Cleavable sites</p> <p data-bbox="1386 341 1512 454">Universal handle</p> <p data-bbox="1554 487 1848 584">a second subsequence (second decoder probe not pictured)</p> <p data-bbox="1386 454 1512 682">Readout Domain 68-80 nt</p> <p data-bbox="1533 682 1680 747">a first subsequence</p> <p data-bbox="1260 714 1428 779">Site-specific linker</p> <p data-bbox="1008 779 1155 828">Antibody</p> <p data-bbox="1029 909 1134 958">Protein</p> <p data-bbox="1428 1088 1890 1120">He-Beechem, Fig. S13 (annotated)</p> <p data-bbox="882 1153 1890 1299">“Protein Detection Probe. The protein detection probe is readily adapted from the RNA detection probe shown in Figure 1. For protein detection, the readout domain contains an additional universal handle and is conjugated to the antibody via a site-specific linker.”</p> <p data-bbox="1428 1347 1890 1380">He-Beechem at 53, notes at Fig. S13</p>

Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		<p>“Following hybridization of ISH probes, slides are washed, assembled into a flow cell, and placed within a fluidic manifold on the SMI instrument for RNA readout and morphological imaging. In RNA assay readout, the tissue is hybridized with 16 sets of fluorescent reporters sequentially; each reporter set contains four single-color reporter pools. The reporters specifically bind to ISH probes during the 16 rounds of reporter hybridization according to the barcode assigned to each gene (Figure S1). After the incubation of each set of reporters, high-resolution Z-stacked images are acquired for downstream analysis.”</p> <p style="text-align: right;">He-Beechem at 5 (emphasis added)</p>  <p style="text-align: right;">He-Beechem, Fig. 1B (annotated)</p>


Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		 <p data-bbox="1459 690 1900 722">He-Beechem, Fig. 1B (annotated)</p> <p data-bbox="882 795 1900 1299">“SMI encoding scheme is designed to assign a unique barcode to each target transcript from a set of 64-bit barcodes (4 color reporters in each readout round over 16 readout rounds) with Hamming distance 4 (HD4) and Hamming weight 4 (HW4). Every barcode is separated by a Hamming distance of at least 4 from all other barcodes to maximally suppress RNA decoding error. Every barcode has a constant HW4, in which each target is ‘on’ in 4 rounds and ‘off’ in 12 rounds during the 16 rounds of reporter hybridization. This ‘on’ and ‘off’ signal barcode design allows for continued expansion to even higher plex since only a subset of RNAs is ‘on’ in any given cycle. Note that the nature of these (on-off) imaging signals represents a ‘deterministic super-resolution imaging system’, and hence each SMI imaging barcode can be located well below the diffraction limit of the imaging system. For each reporter hybridization round, a single reporter can bind to one of the reporter-landing domains on the ISH probe of the gene (Figure 1A and B, Figure S1). The 64-bit encoding scheme with</p>

Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		<p>the HD4 and HW4 yields 1,210 barcodes, from which a subset of 980 barcodes is selected to detect 960 target genes and 20 negative probe controls (Table S1).”</p> <p>He-Beechem at 5 (emphasis added)</p> <p>“The combination of SMI chemistry, hardware, and software enables high-plex spatial RNA profiling within single cells with very high sensitivity. For spatial profiling, the data of fluorescent signals in Z-stacked raw images is transformed into decoded RNA transcript information at their registered spatial location in three dimensions (3D) and assigned to single cells that are identified by segmentation (Figure 1C). During the primary image processing, the fluorescence signals in raw images are transformed into digital output comprising of detected ‘spots’ localized in X-, Y-, and Z-dimensions for subsequent RNA decoding. A spot is identified in an image as an isolated fluorescence signal with the intensity much higher than its neighbors. In biological tissues, however, dense fluorescent spots can be detected at a small spatial region, which may limit the accuracy of resolving each spot. To solve this potential issue, we developed custom image analysis algorithms to process 3D multi-channel image stacks obtained in each field of view (FOV). The key objective of this analytical method is to reduce the multi-dimensional image stack to a single list of individual reporter-binding events. This process is performed across all FOVs, concurrently with image acquisition during cyclic reporter readout. All spots pertaining to a given FOV are collated into a single list of XYZ locations of all individual reporter-binding events. This list is used in the next step to decode the gene-specific barcodes formed by these reporter-binding events.”</p> <p>He-Beechem at 8-9 (emphasis added)</p>
1[c]	(c) comparing said set of signal signatures against set of signal signatures assigned to different	In NanoString’s CosMx SMI workflow, including when performed by NanoString as part of its CosMx Technology Access Program, the set of signal signatures are compared to a list of signal signatures for various analytes include the target analyte to identify the probe and thereby identify the analyte in the sample.

Claim No.	Claim Language	NanoString Spatial Molecular Imaging
	analytes including said analyte, to identify said probe, thereby identifying said analyte immobilized in the sample.	<p>“Each gene is assigned with a unique 16-digit barcode with 4 ‘on’ spots (labeled as ‘1’) and 12 ‘off’ spots (labeled as ‘0’). Each digit of the barcode indicates the presence of the reporter that is associated with the target in the specific reporter hybridization round. ‘1’ means that there is a reporter hybridizing to the ISH probe of the target in that hybridization round, and its color indicates the fluorophore of the hybridized reporter. ‘0’ means that there is no reporter that binds to the target ISH probe in that hybridization round, and the target should be silenced or blank in that round of imaging. For each gene, 4 reporters will sequentially bind to the 4 designated reporter landing domains of the ISH probe throughout the 16 rounds of cyclic reporter readout.”</p> <p style="text-align: right;">He-Beechem at 43, notes at Fig. S1</p>  <p style="text-align: right;">He-Beechem, Figure S1</p> <p>“Analysis workflow to transform raw images to decoded RNA transcripts at subcellular resolution. The workflow includes: 1) three-dimensional primary image processing to identify and register reporter spots, 2) decoding of reporter spots to RNA transcripts with registered X, Y, Z spatial location, 3) outlining of nuclei and cell boundaries with DAPI and antibodies after cyclic</p>

reporter readout for morphology-based cell segmentation, and 4) assigning RNA transcripts to single cells.”

He-Beechem at 8, at Fig. 1C (emphasis added)




He-Beechem, Fig. 1C

“As an example of the final output of the analysis workflow, three identified genes (COL1A1 [yellow], IGHG1 [cyan], and KRT19 [orange]) in FFPE human lung tissue were overlayed with the segmented cells based on their registered spatial information.”

He-Beechem at 8, notes at Fig. 2C

See also He-Beechem, Fig. 4 and notes

Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		<p>“The SMI encoded detection chemistry utilized for RNA can be applied to protein detection by conjugating the oligonucleotide-readout sequences to antibodies (Figure S13). This simple expansion of SMI chemistry for antibody-based detection is enabled by the extremely small size of the oligonucleotides required for SMI 64-bit encoding; only 60 to 80 nucleotides are required to encode all of the multiplexing.”</p> <p>He-Beechem at 19 (emphasis added)</p>



He-Beechem, Fig. 6

Claim No.	Claim Language	NanoString Spatial Molecular Imaging
		<p>“Spatial Subcellular Protein Analysis on SMI. A. Multi-channel overlay of six protein targets detected in a breast cancer biopsy (HER-2 positive invasive carcinoma) from a 108-plex assay (104 encoded targets and 4 morphological markers). B. Enlargement of boxed region in (A). Each decoded marker is visualized along with DAPI (blue) and the morphological markers HER2 (yellow), B7-H3 (light green), CD45 (cyan), CD4 (light purple), PanCK (red), and SMA (white). Scale bars: 50 μm. C. Single-cell expression profiles across 104-encoded protein targets in the sample shown in (A). Colored according to the log2 transformation of the sum of pixel intensities in each cell with a threshold for visualization the geometric mean of the negative controls plus 1.5 standard deviations.”</p> <p>He-Beechem at 21-22, notes at Fig. 6</p>